

Shuanghuanglian Mixture (*Shuanghuanglian Koufuye*)

双黄连口服液

Sample source

Commercially available Shuanghuanglian Mixture

Chemical reference substances

Chlorogenic acid (National Institute for the Control of Pharmaceutical and Biological Products, Batch number: 110753-200212)

Preparation of test solution


Accurately measure 2 mL of the mixture in a 50 mL volumetric flask, add water to volume and mix well.

Chromatographic conditions

- Column: ZORBAX SB C18 4.6×150 mm, 5 µm (883975-902)
- Column temperature: 35 °C
- Mobile phase: methanol-water-glacial acetic acid (20:80:1)
- Detector wavelength: 324 nm
- Flow rate: 1.0 mL/min

Chromatographic system

- Agilent 1200 Series quaternary pump with vacuum degasser
- Agilent 1200 Series high-performance autosampler
- Agilent 1200 Series thermostated column compartment
- Agilent 1200 Series variable wavelength detector
- System control through Agilent ChemStation revision B.01.01



Components	k'	Ret Time (min)	Height (mAU)	Area (mAU*s)	n	USP T _r
Chlorogenic acid	5.433	7.732	76.67	1009.2	8365	0.97