



FIGURE 3.7

With its unique combination of speed and sensitivity, ScanRDI is crucial in microbial detection. (From Eagle Analytical Services. With permission.)



FIGURE 3.8

Dr. Zolner is reviewing lab results after completion of testing at Eagle Analytical Services. (From Eagle Analytical Services. With permission.)

The system is one thing, but the use of the system determines the quality of the reporting of the test results. Eagle Analytical Services performs a three-step procedure on all measurement runs to ensure the quality of its testing results. These steps are as follows: (1) a basic instrumentation performance and calibration test, (2) a zero control test, and (3) a positive control test. Eagle Analytical Services does not issue any reports until all three tests have been successfully run (Figure 3.8).¹⁶

PCCA and the AIB Industry Audit Solution for Compounding Pharmacies

PCCA, inclusive of its sister company Eagle Analytical Services, has a solid base to develop drug-safe audits for compounding pharmacies. These