Where drying in vacuum over a desiccant is directed, a vacuum desiccator, a vacuum drying pistol, or other suitable vacuum drying apparatus shall be used.

6.40.10. Ignite to Constant Weight

"Ignite to constant weight" means that ignition shall be continued at $800 \pm 25^{\circ}$, unless otherwise indicated, until two consecutive weighings, the second of which is taken after an additional period appropriate to the nature and quantity of the residue, do not differ by more than 0.50 mg per g of substance taken.

6.40.20. Dried to Constant Weight

"Dried to constant weight" means that drying shall be continued until two consecutive weighings, the second of which is taken after an additional drying period appropriate to the nature and quantity of the residue, do not differ by more than 0.50 mg per g of substance taken.

6.50. Preparation of Solutions

6.50.10. Filtration

Where a procedure gives direction to "filter" without further qualification, the liquid shall be passed through suitable filter paper or equivalent device until the filtrate is clear. Due to the possibility of filter effects, the initial volumes of a filtrate may be discarded.

6.50.20. Solutions

Unless otherwise specified, all solutions shall be prepared with Purified Water. Solutions for quantitative measures shall be prepared using accurately weighed or accurately meas-

ured analytes (see section 8.20 About).

An expression such as "(1 in 10)" means that 1 part by volume of a liquid shall be diluted with, or 1 part by weight of a solid shall be dissolved in, a sufficient quantity of the diluent or solvent to make the volume of the finished solution 10 parts by volume. For example, a 1 in 10 solution is prepared by diluting 1 mL of a liquid or dissolving 1 g of a solid in sufficient solvent to make 10 mL of the solution.

Auser An expression such as "(20:5:2)" means that the respective numbers of parts, by volume, of the designated liquids shall be mixed, unless otherwise indicated.

6.50.20.1. Adjustments to Solutions

When a specified concentration is called for in a procedure, a solution of other normality or molarity may be used, provided that allowance is made for the difference in concentration and that the change does not increase the error of measurement.

Proportionately larger or smaller quantities than the specified weights and volumes of assay or test substances and Reference Standards may be taken, provided the measure-

ment is made with at least equivalent accuracy.

Unless otherwise indicated, analyte concentrations shall be prepared to within ten percent (10%) of the indicated value. In the case in which a procedure is adapted to the working range of an instrument, solution concentrations may differ from the indicated value by more than ten percent (10%), with appropriate changes in associated calculations. Any changes shall fall within the validated range of the instrument.

When adjustment of pH is indicated with either an acid or base and the concentration is not indicated, appropriate concentrations of that acid or base may be used.

6.50.20.2. Test Solutions

Information on Test Solutions (TS) is provided in the *Test Solutions* portion of the *Reagents, Indicators, and Solutions* section of the *USP-NF*. Use of an alternative Test Solution or a change in the Test Solution used may require validation.

6.50.20.3. Indicator Solutions

Where a procedure specifies the use of an indicator TS, approximately 0.2 mL, or 3 drops, of the solution shall be added unless otherwise directed.

6.60. Units Necessary to Complete a Test

Unless otherwise specified, a sufficient number of units to ensure a suitable analytical result shall be taken.

6.60.10. Tablets

Where the procedure of a Tablet monograph directs to weigh and finely powder not fewer than a given number of Tablets, a counted number of Tablets shall be weighed and reduced to a powder. The portion of the powdered Tablets taken shall be representative of the whole Tablets and shall, in turn, be weighed accurately.

6.60.20. Capsules

Where the procedure of a Capsule monograph gives direction to remove, as completely as possible, the contents of not fewer than a given number of the Capsules, a counted number of Capsules shall be carefully opened and the contents quantitatively removed, combined, mixed, and weighed accurately. The portion of mixed Capsules contents taken shall be representative of the contents of the Capsules and shall, in turn, be weighed accurately.

6.70. Reagents

The proper conduct of the compendial procedures and the reliability of the results depend, in part, upon the quality of the reagents used in the performance of the procedures. Unless otherwise specified, reagents conforming to the specifications set forth in the current edition of *Reagent Chemicals* published by the American Chemical Society (ACS) shall be used. Where such ACS reagent specifications are not available or where the required purity differs, compendial specifications for reagents of acceptable quality are provided (see the *Reagents, Indicators, and Solutions* section of the *USP-NF*). Reagents not covered by any of these specifications should be of a grade suitable to the proper performance of the method of assay or test involved.

Listing of these reagents, including the indicators and solutions employed as reagents, in no way implies that they have therapeutic utility; furthermore, any reference to *USP* or *NF* in their labeling shall include also the term "reagent" or "reagent grade." USP may supply reagents if they otherwise may not be constally commercially available.

wise may not be generally commercially available.

6.80. Equipment

Unless otherwise specified, a specification for a definite size or type of container or apparatus in a procedure is given solely as a recommendation. Other dimensions or types may be used if they are suitable for the intended use.

6.80.10. Apparatus for Measurement

Where volumetric flasks or other exact measuring, weighing, or sorting devices are specified, this or other equipment of at least equivalent accuracy shall be employed.

6.80.10.1. Pipet/Pipette

Where a pipet/pipette is specified, a suitable buret may be substituted. Where a "to contain" pipet/pipette is specified, a suitable volumetric flask may be substituted.

6.80.10.2. Light Protection

Where low-actinic or light-resistant containers are specified, either containers specially treated to protect contents from light or clear containers that have been rendered opaque by application of a suitable coating or wrapping may be used.

6.80.20. Instrumental Apparatus

An instrument may be substituted for the specified instrument if the substitute uses the same fundamental principles of operation and is of equivalent or greater sensitivity and accuracy. These characteristics shall be qualified as appropriate. Where a particular brand or source of a material, instrument, or piece of equipment, or the name and address of a manufacturer or distributor, is mentioned (ordinarily in a footnote), this identification is furnished solely for informational purposes as a matter of convenience, without implication of approval, endorsement, or certification.

6.80.20.1. Chromatographic Tubes and Columns
The term "diameter" refers to internal diameter (ID).

6.80.20.2. Tubing

The term "diameter" refers to outside diameter (OD).