

**Calcium Phosphate Tablets for Cats and Dogs (Direct Compression)**

Bill of Materials			
Scale (mg/tablet)	Item	Material Name	Quantity/1000 Tablets (g)
400.00	1	Dicalcium phosphate	400.00
100.00	2	Wheaten flour	100.00
1.00	3	Citric acid crystalline	1.00
272.00	4	Lactose monohydrate	272.00
QS	5	Flavors	QS
20.00	6	Kollidon® 90F	20.00
4.00	7	Magnesium stearate	4.00

**Manufacturing Directions**

1. Pass all components through an 0.8-mm sieve, mix.

2. Press with medium- to high-compression force (20 kN).
3. Compress into 800-mg tablets, using 12-mm biplanar punches.

**Captopril Tablets (25 mg), Capoten**

Capoten is available in potencies of 12.5, 25, 50, and 100 mg as scored tablets for oral administration. Inactive ingredi-

ents include microcrystalline cellulose, cornstarch, lactose, and stearic acid.

**Captopril Tablets**

Bill of Materials			
Scale (mg/tablet)	Item	Material Name	Quantity/1000 Tablets (g)
25.00	1	Captopril	25.00
91.00	2	Ludipress	91.00
2.00	3	Kollidon CL	2.00
2.00	4	Magnesium stearate	2.00

**Manufacturing Directions**

1. Mix all components, pass through a 0.8-mm sieve, and press with medium-compression force to meet the following specifications.

2. Compress into 122-mg tablets, using 8-mm biplanar punches.