

BATHING-GEL

RECIPE:	% By Weight
A.	
GENAPOL LRO Paste	48.00
HOSTAPUR SAS 60	24.00
GENAPOL AMS	7.00
Oleic acid diethanolamide	1.00
Common salt	2.50
Water	14.50
B.	
Dyestuff	q.s.
Preserving agent	q.s.
Perfume	3.00

## Procedure:

I Melt A at 70C under careful stirring.

II One after another the components of B are added to I at 50C.

Formulation No. A I/5013

BUBBLE-BATH

Clear, middle viscosity

12% active detergent

RECIPE:	% By Weight
A.	
Cocamide DEA	3.00
B.	
Water	15.00
C.	
HOSTAPUR SAS 60	6.00
GENAPOL ZRO Liquid	30.00
Perfume	1.20
Water	42.00
Preserving agent	q.s.
Dyestuff	q.s.
D.	
Common salt	2.80

If GENAPOL ZRP Paste is being used instead of GENAPOL ZRO Liquid, 0.4 times the quantity of GENAPOL ZRO Liquid is diluted with water to the required amount.

## Procedure:

I A is dissolved warm in B.

II One after another the components of C are added to I.

III Finally the viscosity is adjusted with D.

Formulation No. A I/1064

SOURCE: Hoechst Celanese Corp.: Suggested Formulations